



## Immozyme™ hydrolase Kit 2012

ChiralVision's Immozyme hydrolase kit contains 43 immobilized hydrolytic enzymes, ideal for resolutions under aqueous, biphasic and water-free conditions:

- 18 immobilized lipases
- 22 immobilized proteases
- 3 immobilized phospholipases
- 1 gram each
- covalently bound on polyacrylic beads (T2, 150, 250 or 350 µm)
- instructions

*All Immozymes are produced by ChiralVision and available in bulk.*

**Product nr. IMMHYD-4300**

**Pricing kit:** 1875 euro / each\*

\* Shipping included. For ordering information visit the website. Contact us for a US dollar quote.

	<b>Immozyme</b>	<b>source</b>	<b>synonym</b>	<b>class</b>
1	IMMCALA-T2-350	<i>Candida antarctica, A</i>	CaLA <sup>1</sup>	lipase
2	IMMCALB-T2-350	<i>Candida antarctica, B</i>	CaLB <sup>1</sup>	lipase
3	IMMCALBY-T2-350	<i>Candida antarctica, B</i>	CaLBY <sup>2</sup>	lipase
4	IMMRML-T2-350	<i>Rhizomucor miehei</i>	Lipozyme <sup>1</sup>	lipase
5	IMMTLL-T2-350	<i>Thermomyces lanuginosa</i>	Lipolase <sup>1</sup>	lipase
6	IMMCRL-T2-350	<i>Candida rugosa</i>	CRL <sup>3</sup>	lipase
7	IMMABC-T2-350	<i>Pseudomonas cepacia</i>	Amano PS <sup>4</sup>	lipase
8	IMMAPF-T2-350	<i>Pseudomonas fluorescens</i>	Amano AK <sup>4</sup>	lipase
9	IMMARO-T2-350	<i>Rhizopus oryzae</i>	Amano F-AP15 <sup>4</sup>	lipase
10	IMMAMJ-T2-150	<i>Mucor javanicus</i>	Amano M <sup>4</sup>	lipase
11	IMMANA-T2-150	<i>Aspergillus niger</i>	Lipase A <sup>4</sup>	lipase
12	IMMRNA-T2-150	<i>Rhizopus niveus</i>	Lipase <sup>4</sup>	lipase
13	IMMASMQ-T2-150	<i>Alcaligenes sp.</i>	lipase QLM <sup>6</sup>	lipase
14	IMMRES-T2-150	Resinase HT	Resinase HT <sup>1</sup>	lipase
15	IMMLIPX-T2-150	Lipex 100L	Lipex 100L <sup>1</sup>	lipase
16	IMML51-T2-350	Novozymes 51032	51032 <sup>1</sup>	lipase
17	IMMCCMO-T2-150	<i>C. cylindracea (C. rugosa)</i>	lipase OF <sup>6</sup>	lipase
18	IMMAULI-T2-350	<i>Bacillus subtilis</i>	lipase <sup>7</sup>	lipase
19	IMMALC-T2-250	<i>Subtilisin from Bacillus sp.</i>	Alcalase <sup>1</sup>	protease
20	IMMSAV-T2-250	<i>Subtilisin from Bacillus sp.</i>	Savinase <sup>1</sup>	protease

	<b>Immozyme</b>	<b>source</b>	<b>synonym</b>	<b>class</b>
21	IMMEVE-T2-250	<i>Subtilisin from Bacillus sp.</i>	Everlase <sup>1</sup>	protease
22	IMMESP-T2-250	<i>Subtilisin from Bacillus sp.</i>	Esperase <sup>1</sup>	protease
23	IMMCAR-T2-150	<i>Mucor miehei</i>	Acid protease <sup>1</sup>	protease
24	IMMP6-T2-250	<i>Bacillus licheniformis</i>	Protex 6L <sup>8</sup>	protease
25	IMMP7-T2-150	<i>Bacillus amyloliquefaciens</i>	Protex 7L <sup>8</sup>	protease
26	IMMP8-T2-250	<i>Bacillus licheniformis</i>	Protex 8L <sup>8</sup>	protease
27	IMMP14-T2-150	<i>Geobacillus sp. (Thermolysin)</i>	Protex 14L <sup>8</sup>	protease
28	IMMP15-T2-150	<i>Trichoderma reesei</i>	Protex 15L <sup>8</sup>	protease
29	IMMP26-T2-150	<i>Aspergillus niger</i>	Protex 26L <sup>8</sup>	protease
30	IMMP30-T2-150	<i>Bacillus subtilis</i>	Protex 30L <sup>8</sup>	protease
31	IMMP40L-T2-150	<i>Bacillus subtilis</i>	Protex 40L <sup>8</sup>	protease
32	IMMP50-T2-150	<i>Aspergillus oryzae var.</i>	Protex 50 FP <sup>8</sup>	protease
33	IMMP51-T2-150	<i>Aspergillus oryzae</i>	Protex 51 FP <sup>8</sup>	protease
34	IMMP89-T2-250	<i>Bacillus subtilis</i>	Protex 89L <sup>8</sup>	protease
35	IMMAUAC-T2-150	<i>Aspergillus niger</i>	Acid protease <sup>7</sup>	protease
36	IMMAUNE-T2-150	<i>Bacillus subtilis</i>	Neutral protease <sup>7</sup>	protease
37	IMMAUAL-T2-250	<i>Bacillus subtilis</i>	Alkaline protease <sup>7</sup>	protease
38	IMMPAP-T2-150	<i>Carica papaya</i>	Papain <sup>3</sup>	protease
39	IMMBRO-T2-150	<i>pineapple stem</i>	Bromelain <sup>3</sup>	protease
40	IMMFIC-T2-150	<i>fig tree latex</i>	Ficin <sup>3</sup>	protease
41	IMMLMPL-T2-150	LysoMax	phospholipase A2 <sup>8</sup>	phospholipase
42	IMMGZPL-T2-150	G-Zyme G999	lyso-phospholip. <sup>8</sup>	phospholipase
43	IMMLEPL-T2-150	<i>Lecitase ultra</i>	phospholipase <sup>1</sup>	phospholipase

Registered trademark from: 1) Novozymes; 2) c-LEcta 3) Sigma-Aldrich 4) Amano; 5) Fluka; 6) Meito 7) Aum Enzymes, 8) Genencor

Orders can be send to: [sales@chiralvision.com](mailto:sales@chiralvision.com)

Info: <http://chiralvision.com/order.htm>